

SHORT PLANT

Key Features

Growth height: **20**
inches, 510 mm

Growth area: **11.3**
ft² (1.05 m²)

Light intensity: **550**
μmols/m²/sec

Lighting: **LED**
Light emitting diode

Airflow: 
Horizontal

***Max loading:** **≈48**
Pots, 6" (150 mm) diameter

Horizontal airflow over multiple shelves optimizes growth area for shorter plants, like Arabidopsis. Ideal for research in propagation, genetics, physiology and other moderate light experiments.

