

MTPS FLEX 2-TIER

Key Features

Growth height: **36**
inches, 914 mm

Growth area: **48-288**
ft² (4.5-26.8 m²)

Light intensity: **500**
 $\mu\text{mol}/\text{m}^2/\text{sec}$

Lighting: **FL & H/IN**
Fluorescent &
Halogen/Incandescent

Airflow: 
Horizontal

*Max loading: **≈192-6,912**
Pots, 6" (≈150 mm)
diameter

MULTI-TIER WALK IN

The MTPS Flex allows the room to be converted from a tall plant/high light/single tier configuration to a short plant/low light/two or three-tier configuration. The flexible design means there is no need for special trades or services to adjust shelf height and airflow as research needs change. Available in 72 ft² (6.7 m²) increments.

