

MTR30

PLANT GROWTH CHAMBER

Key Features

Growth height: **25**
inches, 635 mm

Growth area: **30**
ft² (2.79 m²)

Light intensity: **550**
μmols/m²/sec

Lighting: **FL & H/IN**
Fluorescent &
Halogen/ Incandescent

Airflow: 
Upward

***Max loading:** **≈360**
Seed cells
≈2" x 2" (5x5cm)

The MTR30 contains two shelves and is suitable for a range of applications including the growth of short to medium height plants that require moderate light intensity.

