

PGR15

PLANT GROWTH CHAMBER

Key Features

Growth height: **57**
inches, 1450 mm

Growth area: **16**
ft² (1.5 m²)

Light intensity: **1,000**
 $\mu\text{mol}/\text{m}^2/\text{sec}$

Lighting: **LED**
Light Emitting Diode

Airflow: 
Vertical

*Max loading: **≈65**
Pots, 6" (150 mm) diameter

Single tier chamber with high light and counterbalanced lighting enables height adjustment between the lamp canopy and the research material.

